

Store Section
Dr. BRA IRCH

Ref No: F. No. 25/IRCH/SSK/L.O/16-17

Date: 17/05/2016

Subject: Purchase of Channel Pipettes & Tubes, Mfg.by eppendorf AG Germany for Lab.Onco. Dept. (Dr. BRAIRCH)

The request letter received from Dr. Anita Chopra duly forwarded by Dr. Rajive Kumar (Prof. & Head Lab Oncology) for the procurement of Channel Pipettes & Tubes, on proprietary basis from M/s eppendorf AG Germany

Sr. No.	Description Item /Cat No#
1.	3122 000.051 Research Plus, 8-channel 30-300µl Fully Autoclavable Spring Loaded Tip Cone: Improved ergonomics, Tight fit to the tip (not in 5ml and 10ml pipettes) Adjustment opening: Adjust your pipettes to a specific liquid and volume. Control Button: Very low operating force, Color indicates pipette volume, Positioned for perfect ergonomics. Removable independent channels: Gain flexibility for individual needs Ejector: Very low operating force, positioned for perfect ergonomics. Volume Display: 4 Digits, Magnifying shape. Eppendorf Perfect Piston System: Ultra light system made of Forton. Quick Connection Clip: Remove lower part easily
2.	3880 001.018 IsoTherm-system Starter Set, for 1.5/2.0 ml tubes Application Transport and slow thawing of deep-frozen samples Product features Holds 24 tubes (Safe-Lock 1.5/2.0 ml) IsoPack (cold pack), white for 0°C and blue for -21°C, cooled to refrigeration temperature in freezer compartment Isothermal times of < 6 h (0°C, white) and < 3 h (-21°C, blue) No more melting ice baths (danger of contamination) All system components can be used separately Racks stackable, autoclavable, markable, floatable and can be centrifuged in microtiter-plate rotor

The user department has submitted proprietary Article certificate in addition to one submitted by company.

The above documents are being uploaded on AIIMS WEBSITE for open information to invite objection, comments, if any, from any manufacturer regarding proprietary nature of the item within 15 days from the date of issue/uploading of the notification giving reference F. No 25/IRCH/SSK/L.O/16-17. The comments should be sent to stores officer, Dr. BRAIRCH AIIMS on or before 03/06/2016, up to 4:30 p.m failing which it will be presumed that any other vendor/company is having no comment to offer and case will be decided purely on merits.


Stores Officer
Dr. BRA IRCH

Encl:

- PAC Certificate (enclosed)

Date: 22nd October, 2012

TO WHOM IT MAY CONCERN

We here by certify that "Eppendorf Research Plus Pipette" is the proprietary in nature because of its specification & Special features like (Spring loaded TIP Core; Adjustment window to handle viscous solutions; Fortron Piston) manufactured by M/s Eppendorf AG Germany.

To the best of our knowledge; no other manufacturer world wide has similar/ identical Specification.

Yours faithfully,

Dr. Andrés Jarrin
Area Sales Manager
Eppendorf AG
Barkhausenweg 1
22339 Hamburg
Germany

eppendorf^{AG}
Global Sales
Barkhausenweg 1
22339 Hamburg
Germany

Phone: +49 40-5 38 01-172
Fax: +49 40 5 38 01-840
jarrin.a@eppendorf.de



डॉ. अनीता चोपड़ा, एम्बी. डीएस/Dr. Anila Chopra, MD, DM
सहायक आचार्य/Assistant Professor
प्रयोगशाला अतुर्विज्ञान एकाद/ Laboratory Oncology Unit
डॉ. भा.रा.अ.सं.सो.कै.अ. अ.भा.आ.सं., नई दिल्ली-29
Dr. BRAIRCH, AIIMS, New Delhi-110029

Vorstand:
Dr. Dirk H. Ehlers (Vorsitzender)
Dipl.-Kfm. Detmar Ammermann
Dr. Heinz Gerhard Kohn
Dr. Michael Schroeder

Vorsitzender des Aufsichtsrats:
Dipl.-Wirtschaftsling. Klaus Fink

Eppendorf AG
Barkhausenweg 1, 22339 Hamburg
Postanschrift: 22331 Hamburg
Tel.: (040) 5 38 01-0
Fax: (040) 5 38 01-556
www.eppendorf.com

Amtsgericht: Hamburg HRB 76249

Deutsche Bank AG
Kto. 0202200 / BLZ 200 700 00
IBAN DE51 2007 0000 0020 2200 00
S.W.I.F.T. DEUT DE 33

HSBC Trinkaus & Burkhardt K&A
Kto. 600547003 / BLZ 300 300 80
IBAN DE56 3003 0800 0800 5470 08
S.W.I.F.T. TUBD DE 33

Hamburger Sparkasse
Kto. 1280118447 / BLZ 200 505 50
IBAN DE32 2005 0550 1280 1184 47
S.W.I.F.T. HASP DE 33