All INDIA INSTITUTE OF MEDICAL SCIENCES ANSARI NAGAR, NEW DELHI – 110029 Department of Surgical Discipline (Stores Section)

Ref. No. 15/Surgery/Proprietary/2019-20/St.

Dated: 05.12.2019

Subject: Purchase of Advance Table Mounted Abdominal Retractor System Surgical Exposure System on proprietary basis rates – reg.

The request has been received from the Surgeon of Department of Surgical Discipline for purchase of subject item from M/s Thomson Surgical Instruments through their authorized distributors M/s Ultimate Medical Devices on proprietary basis. The proposal submitted by M/s Ultimate Medical Devices , 19, 2nd Floor, Guru Nanak Market, Lajpat Nagar-IV, New Delhi PAC certifications are attached herewith.

The above document are being uploaded on AIIMS Website for open information to prospective manufactures to submit their objections/comments, (if any) regarding proprietary nature of the equipment/item within 15 days giving reference No. 15/Surgery/Prop./2019-20/St. The comments should be received in the office of Sr. Stores Officer (Surgery), Department of Surgical Disciplines, Center for Medical Education & Technology (CMET), Educational Hall, 1st Floor, (Near B.B. Dikshit Library), AIIMS, Ansari Nagar, New Delhi-110029 on or before 23/12/2019 up to 04:00 p.m., failing which, it will be presumed that any other vendor is having no comment to offer and case will be decided on merits.

Sr. Stores Officer (Surgery)

Encl: Related documents enclosed.

- 1. PAC Certificate enclosed.
- 2. Specification of equipment.



18th May 2018

PROPRIETORY CERTIFICATE

This is to certify that all Thompson Surgical products, namely the Thompson Retractors with Quick-Frame and Cam-Type Joints, as well as the list below, are the proprietary products of Thompson Surgical and are exclusively manufactured by Thompson Surgical Instruments, Inc., at 10170 E. Cherry Bend Road, Traverse City, Michigan, 49684 – USA and any of the Thompson product cannot be imitated or copied by any other manufacturer/firm, or person anywhere in world.

Crank retractor handle US 8696562 B2 Crank retractor handle US 7338442 B2 Low profile surgical retractor US 8753272 B2 Adjustable rail clamp with clamp locking device US 8617064 B2 Surgical retractor angling device US 8360971 B2 Surgical retractor with locking blade US 8257255 B2 Surgical retractor having locking interchangeable blades US 5984865 A Adjustable rail clamp US 8100827 B2 Posterior lumbar retractor system US 9241619 B2 Chirurgisches Retraktorsystem DE 19936570 B.1 Cross bar for a surgical retractor system US 6511423 B2 Retractor system for anterior cervical spine surgery US 9320506 B2 Universal scissors joint apparatus US 7749163 B2 Angling surgical retractor apparatus and method of retracting anatomy US 5902233 A Cam tightened universal joint clamp US 6017008 A Retractor System having an articulated frame WO2016191683A1 Surgical Retractor having clamping mechanism US20170042527A1 Surgical retractor with asymmetric blade US20170065268A1

The entire range of Thompson products is promoted by Últimate Medical Devices having its office at 19, 2nd Floor, Guru Nanak Market, Lajpat Nagar – IV, New Delhi – 110 024, exclusively in North & North-East India.

For and on behalf of

Thompson Surgical Instruments, Inc,

Jennifer Dixon Office Manager

International Customer Service

Thompson Surgical Instruments, Inc.





<u>Specifications for Advance Multi-Purpose Table Mounted</u> <u>Exposure System</u>

- Advanced multipurpose table mounted retractor system for Upper abdominal, Hepatobiliary, colo-rectal and pelvic surgery.
- Should have bilateral frame to provide exceptional costal margin retraction.
- Should have various detachable blade options in different configurations for retraction of abdominal wall, solid &hollo viscera. Blades like Balfour, Richardson, Malleable, Harrington. Deaver, Kelly etc. along with pediatric blades.
- All blades should be detachable from blade holder & simultaneously different blade should be interchangeable to any/ all blade holders.
- Blade holder should have inbuilt angling system on the blade holder.
- Blade & blade holder should have fixed retraction feature to hold organ
- Blade holder rods should have ratchet mechanism in built for retraction.
- Should be of high quality steel & US FDA approved
- Physical demonstration must.
- Two years warranty.

| SL | DESCRIPTION | | |
|----|--|-----|--|
| 1 | ELITE II RAIL CLAMP JOINTS 18" | 2 | |
| 2 | CROSSBAR WITH 2 CAM JOINTS - WIDE | | |
| 3 | BILATERAL CROSSBAR HINGED 34" - 111/2x81/2x111/2 | | |
| 4 | BILATERAL CROSSBAR | | |
| 5 | ANGLED ARM 24" - 8X16 | 1 2 | |
| 6 | ANGLED ARM 18" - 8X10 | 2 | |
| 7 | ANGLING BLADE HOLDER 8" | 2 | |
| 8 | ANGLING MICRO ADJUSTABLE BLADE HOLDER 73/4" | 2 | |
| 9 | ANGLING MICRO ADJUSTABLE BLADE HOLDER 10" | 4 | |
| 10 | ADAPTER FOR HAND HELD | 1 | |
| 11 | ANGLING WRENCH | 1 | |
| 12 | CONTAINER WITH SILICON MAT | 2 | |
| 13 | BALFOUR BLADE WITHOUT LIPS SIZE - 70MM X 76MM | 2 | |
| 14 | BALFOUR BLADE WITHOUT LIPS SIZE - 83MM X 76MM | 2 | |
| 15 | BALFOUR BLADE SIZE - 83MM X 57MM | 1 | |
| 16 | KELLY BLADE SIZE - 64MM X 76MM | 1 | |
| 17 | KELLY BLADE SIZE - 76MM X 89MM | 1 | |
| 18 | SOFT MALLEABLE BLADE SIZE - 25MM X 102MM | 1 | |
| 19 | MALLEABLE BLADE SIZE - 51MM X 203MM | 1 | |

| 20 | ST. MARKS BLADE SIZE - 64MM X 178MM | 1 |
|----|--|---|
| 21 | HARRINGTON BLADE SIZE - 64MM X 152MM | 1 |
| 22 | DEAVER BLADE SIZE - 51MM X 127MM | 1 |
| 23 | RICHARDSON BLADE SIZE - 51MM X 127MM | 1 |
| 24 | RADIO RICHARDSON BLADE SIZE - 51MM X 178MM | 1 |
| 25 | MALLEABLE FINGER BLADE SIZE - 6" | 1 |

Technical Specification for Thompson Retractor HOD , Department of Surgery

| SL | ARTICLE | DESCRIPTION | QTY |
|----|-----------|---|-----|
| | | * / | |
| 1 | 41903AC | Elite II Rail Clamp w/2 Cam Joints 18" | 2 |
| 2 | 41900BCW | Crossbar w/2 Cam Joints - Wide | 1 |
| 3 | 41900PH | 34" Bilateral Crossbar Hinged 11½x8½x11½ | 1 |
| 4 | 41900F | Bilateral Crossbar | 1 |
| 5 | 44124N | 24" Angled Arm 8" x 16" | 2 |
| 6 | 44118N | 18" Angled Arm 8" x 10" @45° | 2 |
| 7 | SL42126WA | S-Lock Handle Angling 8" | 2 |
| 8 | SL45007CA | S-Lock Micro-Adj II Clip-on Angling 73/4" | 2 |
| 9 | SL45006CA | S-Lock Micro-Adj II Clip-On Angling 10" | 4 |
| 10 | SL42128F | S-Lock Handle Held Adapter | 1 |
| 11 | 43032 | Angling Head Wrench | 1 |
| 12 | 50000BL | Instrument Case 26 x 10 x 5" | 2 |
| 13 | SO46140C | Balfour w/o Lips 70mm x76mm(2¾" x 3")-SL | 2 |
| 14 | S046141C | Balfour w/o Lips 83mm x76mm (3¼" x3")-SL | 2 |
| 15 | S046241A | Balfour 83mm x 57mm (31/4" x 21/4")-SL | 1 |
| 16 | S046129 | Kelly 64mm x 76mm (2½" x 3")-SL | 1 |
| 17 | S046130 | Kelly 76mm x 89mm (3" x 3½")-SL | 1 |
| 18 | SL46118P | Soft Malleable 25mm x 102mm (1" x 4")-SL | 1 |
| 19 | SL46122 | Malleable 51mm x 203mm (2" x 8") - SL | 1 |
| 20 | S046188 | St Marks 64mm x 178mm (2½"x7") w/Lip-SL | 1 |
| 21 | SL46160 | Harrington 64mm x 152mm (2½" x 6")-SL | 1 |
| 22 | SL46174 | Deaver 51mm x 127mm (2" x 5") - SL | 1 |
| 23 | S046151 | Richardson 51mm x 127mm (2" x 5")-SL | 1 |
| 24 | S046153ET | Radio Richardson 51mm x 178mm - SL | 1 |
| 25 | SL46191 | Malleable Finger 6" - SL | 1 |